





15ml Centrifuge tube with Cap



Available for any kind of operation and sample storage.
Spins to 10000xg
Autoclavable (121°C / 15min)

Item :
Natural : SG6555-000
Amber : SG6555-900

Package :
25 tubes / bag , 20 bags / cs ,(500 tubes / cs)

Made by Special materials (Polypropylene)
Occurs writing panel and easy-to-read graduations(1 ml per line)

Only for research.
DNase & RNase Free, Pyrogen Free, Endotoxin Free, Human DNA Free.



Cat No	specification
6555-000	Natural 15ml Centrifuge tube with Cap
6555-900	Amber 15ml Centrifuge tube with Cap