




Smart Pipet



Smart pipet is an electronic pipette controller designed for use with graduated and volumetric glass and plastic pipets. The ergonomic design make it very easy handle and operate. The pipetting speed can be controlled by push buttons or continuous speed control regulator. The dispense can be driven by gravity or blow-out with push button.

The re-chargeable Lithium battery is the most advanced battery technology available in the market. The pipette controller holder with optional magnet enables the users to hang the unit on the wall or laminar flow. UV-resistant plastic materials enables the unit to endure short period of time of UV light.

- * **Battery Lidicator**
- * **Helix Switch Close-knit**
- * **UV-resistant housing**
- * **Lithium Battery**
- * **Blow-out/gravity Dispense switch**
- * **Continuous Speed Control Wheel**

Cat No	specification
1001-000	Smart Pipet
1002-000	Hydrophobic filter
	10` s/bag

- ◆ Powerful Lithium battery, No memory effect
- ◆ Continuous speed control wheel
- ◆ Blow-out and gravity dispense selection
- ◆ Dangerous fume or liquid protective hole
- ◆ Single-handed Operation , Self-inverted on the bench
- ◆ UV-resistant housing , Complies with CE and RoHS
- ◆ Light-weighted; 167g, including a battery
- ◆ Other bands to serve as reference points.