

Smart Pellet pestles Cordless motor



- ◆ Hand-operated or motor-driven grinders for resuspending pellets or disrupting soft tissue in microcentrifuge tubes.
- ◆ Its high-speed vortexing action completely homogenizes cell tissue and pellets within seconds.
- ◆ Designed to homogenize cells and plant tissue or to re-suspend protein and DNA pellets.

- ◆ Pestle ends are specially designed to mate with 0.5 mL or 1.5 mL microcentrifuge tubes.
- ◆ Shaft is convenient length and diameter for gentle manual back-and-forth rotation or more vigorous motorized operation.
- ◆ Optional battery-powered motor (with two AA batteries) can be used with pestles and microcentrifuge tubes.
- ◆ Motor is lightweight and greatly reduces hand fatigue.

Cat No	specification
5140-000	Smart Pellet pestles Cordless motor

