

Autoclavable bag Sterilisation Indicator Tape



Cat No	specification	
5136-500	13*17", Orange, Adhesives, printed "BIOHARZARD"	100's/bag
5130-500	24*30", Orange, Tie wire, printed "BIOHARZARD"	100's/bag
5136-000	33*45cm, Transparent, Adhesives, non-printed	100's/bag
5139-000	60*90cm, Transparent, Tie wire, non-printed	100's/bag
7161-000	Sterilisation Indicator Tape, Size : W12mm x L50M	30Roll/box

Our rich industry experience enables us to present an extensive array of Bio Medical Waste Bag to our clients. Used for collecting and disposing the wastage of pharmaceutical industries and laboratories, this assortment of waste bag is fabricated using superior grade PVC, We have two kind of printed and non-printed logo Bio Medical Waste Bag. Our offered biomedical waste collection bag is highly appreciated in the market for to its durability and perfect finish and is available in various colors and sizes to meet the requirements of clients.

CAUTION : laboratory safety regulations require that all biohazardous materials be autoclaved prior to disposal.

NOTE : Do not seal bag tightly during autoclaving.
DO NOT EXCEED 260°F or 125°C