

Cap Opener & Cap Holding Tabs & Floating Pad

Cap Holding Tabs :

- ◆ These holding tabs secure lids on either 1.5ml, 1.7ml, 2.0ml standard microcentrifuge tubes to provide extra security and facilitate tubing removal.
- ◆ Material : Polypropylene

Cap Opener :

- ◆ Can be use for 1.5ml, 1.7ml, 2.0ml microcentrifuge tubes, Autoclavable.
- ◆ Material : Polypropylene

Floating Pad :

- ◆ This floating pad can be use float in your water bath.

Cat No	specification	
5121-700	Floating Pad	5/ bag
5111-800	Cap Holding Tabs (Assorted colors)	100/ bag
5111-610	Cap Holding Tabs (Pink colors)	100/ bag
5112-000	Cap Opener & Cap Holding Tabs & Floating Pad set	SET/ bag
	Removable or in combination, This package include :	
	Cap Opener × 1 、 Cap Holding Tabs × 8 、 Floating Pad × 1	

