

## Electroporation cuvettes



Cat No	specification	
4301-900	1mm Electroporation cuvettes, Brown	1/bag, 25/pkg
4301-200	2mm Electroporation cuvettes, Green	1/bag, 25/pkg
4301-100	4mm Electroporation cuvettes, Blue	1/bag, 25/pkg

- ◆ For Electroporation of DNA, RNA, & protein in bacteria, yeast, plant & animal cells.
- ◆ Compatible with most major brands of electroporation instrumentation
- ◆ Correct positioning in the cuvette holder guaranteed due to design
- ◆ Plastic cuvette with aluminum electrodes 3 sizes available 100, 400, and 800  $\mu$ l (1 mm, 2 mm, and 4 mm gap width)
- ◆ Individually sealed and sterilized by UV/gamma irradiation
- ◆ Color-coded for easy identification