

## PP Disposable Pellet Pestles



Cat No	specification	
5142-000	PP Disposable Pellet Pestles For 0.5ml tube	100/box
5141-000	PP Disposable Pellet Pestles For 1.5ml tube	100/box
5140-000	Motor (include 0.5ml or 1.5ml sample pestle 100pcs)	SET/box

- ◆ Material : Polypropylene
- ◆ To grind the sample and DNA pellet in 1.5ml or 0.5ml centrifuge tube.
- ◆ Pestle length is 1 cm

