

Pipetting Reservoir



Cat No	specification
5143-000	Pipetting Reservoir , 5` s/bag, 200` s/cs
5143-010	Pipetting Reservoir , 5` s/bag, 20` s/box
5144-000	Pipetting Reservoir Sterile · Single Pakg. 100` s/cs
5144-010	Pipetting Reservoir Sterile · Single Pkg. , 5` s/bag, 200` s/cs

Reagent Reservoir, autodavable reservoir , filling volume 50-60 ml, Convenient autodavable reservoir. Easy-to-use loading of multichannel pipettes. The Reagent Reservoir is an easy way of collecting fluids distributed by a multichannel pipette.

The Reagent Reservoir is kept clean by means of an autoclave; the reservoir is designed to resist high pressured steam at temperatures up to 121 degrees C during sterilization in an autoclave.