

## Square Sealing Mats



Square sealing mats evaporation in a wide range of temperatures in PCR and storage applications. Square sealing mats are chemically resistant, reusable sealing mats, that provide a high degree of confidence. For storage applications involving more aggressive solvents, they contain no adhesives or extractable to contaminate samples. Made of research-grade silicone, they can be autoclaved, they can also be reused.

Contains no adhesives, no extractables to contaminate samples

- ◆ Tight seal to the plate, no sample evaporation or well-to-well contamination
- ◆ Usable within wide temperature range from -80°C up to +121°C
- ◆ High contrast alphanumeric grid for easy well identification
- ◆ For use with 96 square well deep well plates

Cat No	specification
7152-000	Square Sealing Mats

5's/bag, 50/case