

96Well Foil Seal & Plate Seal (Adhesive/Thermal Seal)



Foil Adhesive Seal



Adhesive Plate Seal



Thermal Seal



Roll Wheel

Cat No	specification	
7112-000	96Well Adhesive Plate Seal (Piece)	100/bag
7131-000	96Well Adhesive Foil Seal (Piece)	100/bag
7141-000	Thermal Seal Foil , 80mm*250M (Roll)	250M/roll
7141-000	Thermal Seal Foil , 80 x 135 mm (Piece)	100/bag
	Roll Wheel	

- ◆ Strong adhesion prevents evaporation during thermal cycling or sample storage
- ◆ Sample is accessed by piercing foil with a pipette tip
- ◆ No adhesive build-up on tip after piercing
- ◆ End tabs is easy for seal and removal application
- ◆ The Roller tools may bind the plate easier and closer (Product NO.6241)