

QPCR 8-strip tube & caps / individually strip tube



4032-000



Cat No	specification	
4031-000	Optically Clean Strips tube & cap	120` s/box , 10 box/ cases
4231-000	Optically Clean Strips cap	250` s/box , 10 box/ cases
4032-000	8-Strip PCR Tube with Individually Attached Optically Clear Flat Cap	120` s/box , 10 box/ cases

Each Ultra Clear Cap has a highly polished flat surface, providing an optically clear window. The caps are mounded from a specially formulated polypropylene that enhances signal detection by minimizing autofluorescence and maximizing transmission.

- ◆ SG4031-000 type : 8-cap strip and tube
- ◆ SG4032-000 type : 8-Strip PCR Tube (Natural) with Individually Attached Optically Clear Flat Cap, Tubes attached by One Strand
- ◆ Specially formulated ultra clear
- ◆ Polypropylene
- ◆ Highly polished
- ◆ Flat cap design