

## 96 well Non-Skirt PCR Plate

- ◆ Suitable for 0.2ml thermal cycler blocks
- ◆ 96-well PCR plate
- ◆ 0.3ml maximum well volume (when used with adhesive and heat seals)
- ◆ Cut-away corner
- ◆ Fits all 0.2ml thermal cyclers
- ◆ Clear well bottom for sample visibility
- ◆ Highly rigid plate deck

Cat No	specification	
4061-000	0.2ml 96 well Non-Skirt PCR Plate , Nature	25` s/box , 5box/cs
4061-300	0.2ml 96 well Non-Skirt PCR Plate , Yellow	25` s/box , 5box/cs
4061-200	0.2ml 96 well Non-Skirt PCR Plate , Green	25` s/box , 5box/cs
4061-100	0.2ml 96 well Non-Skirt PCR Plate , Blue	25` s/box , 5box/cs
4061-400	0.2ml 96 well Non-Skirt PCR Plate , Purple	25` s/box , 5box/cs
4061-W00	0.2ml 96 well Non-Skirt PCR Plate , White	25` s/box , 5box/cs
4061-700	0.2ml 96 well Non-Skirt PCR Plate , Black	25` s/box , 5box/cs

