

96 well Semi-Skirt PCR Plate



- ◆ Suitable for 0.2ml thermal cycler blocks
- ◆ 96-well PCR plate
- ◆ 0.3ml maximum well volume (when used with adhesive and heat seals)
- ◆ Cut-away corner
- ◆ Clear well bottom for sample visibility
- ◆ Highly rigid plate deck

Cat No	specification	
4062-000	0.2ml96 well Semi-Skirt PCR Plate , Nature	25 ` s/box , 5box/cs
4062-300	0.2ml96 well Semi-Skirt PCR Plate , Yellow	25 ` s/box , 5box/cs
4062-200	0.2ml96 well Semi-Skirt PCR Plate , Green	25 ` s/box , 5box/cs
4062-100	0.2ml96 well Semi-Skirt PCR Plate , Blue	25 ` s/box , 5box/cs
4062-400	0.2ml96 well Semi-Skirt PCR Plate , Purple	25 ` s/box , 5box/cs
4062-W00	0.2ml96 well Semi-Skirt PCR Plate , White	25 ` s/box , 5box/cs
4062-700	0.2ml96 well Semi-Skirt PCR Plate , Black	25 ` s/box , 5box/cs