


## Combined type Four Faces Rack



### 15ml Four Faces Rack:

- 50ml(  $\phi$  30mm) 2well
- 15ml(  $\phi$  17mm) 16+6well
- 1.5ml(  $\phi$  12mm) 15+9+6well
- 0.5ml(  $\phi$  8mm) 28well



### 50ml Four Faces Rack:

- 50ml(  $\phi$  30mm) 5+2well
- 15ml(  $\phi$  17mm) 6well
- 1.5ml(  $\phi$  12mm) 26+6+8well
- 0.5ml(  $\phi$  8mm) 13well

Cat No	Color	specification	
5032-600	Pink	50ml Combined type Four Faces Rack	100's/cs
5031-600	Pink	15ml Combined type Four Faces Rack	100's/cs
5032-100	Blue	50ml Combined type Four Faces Rack	100's/cs
5031-V00	Purple	15ml Combined type Four Faces Rack	100's/cs
5032-700	Red	50ml Combined type Four Faces Rack	100's/cs
5032-800	Yellow	50ml Combined type Four Faces Rack	100's/cs
5032-900	Green	50ml Combined type Four Faces Rack	100's/cs

### Combined type Faces Rack Description :

- ◆ PP test tube rack can be used for multipurpose in four face different test tubes diameter.
- ◆ Polypropylene racks and boxes suitable for use with 0.5ml · 1.5ml · 1.7ml · 2.0ml · 15ml · 50ml tubes
- ◆ Multi-use depending on volume required.
- ◆ Liquid volume can be easily told.