



# 96-Place double-sided Rack ( For Microcentrifuge tube / Drying Rack )




### Front

Rack for 96-Place  
Microcentrifuge tube  
It can hold 1.5 and 2ml  
tubes and cryo-vials.  
Wide-space of holes.  
Easily take and operate.



### Back

According to different uses,  
place the removable column,  
it can be Drying rack or  
another size of centrifuge rack.  
Depends on the column you  
placed.  
(Columns aren't included.)



Cat No	96-Place double-sided Rack (For Microcentrifuge tube / Drying Rack) 20` s/cs				
<b>2500-000</b>	Nature	96-Place double-sided Rack	<b>2500-000R</b>	Red	96-Place double-sided Rack
<b>2500-000B</b>	Blue	96-Place double-sided Rack	<b>2500-000Y</b>	Yellow	96-Place double-sided Rack
<b>2500-000G</b>	Green	96-Place double-sided Rack	<b>2501-000</b>	White	Columns 100` s/bag