

80 Well Double Sided Freezer Rack

- ◆ 1.5/2.0ml tube, and the other side can be operate the 0.2/0.5 ml tube
- ◆ Can in under the low temperature environment operation
- ◆ Have every kind of color, the convenient laboratory
- ◆ Economical space
- ◆ classification manages
- ◆ Material : Polypropylene

Cat No	specification	
5021-000	80 well double Sided freezer rack (natural)	100/case
5021-300	80 well double Sided freezer rack (yellow)	100/case
5021-400	80 well double Sided freezer rack (purple)	100/case
5021-100	80 well double Sided freezer rack (blue)	100/case
5021-B00	80 well double Sided freezer rack (light blue)	100/case
5021-200	80 well double Sided freezer rack (green)	100/case
5021-500	80 well double Sided freezer rack (orange)	100/case

