





96 Place PCR Tube Rack



Cat No	specification	
5011-400	96 well PCR tube Racks (Purple)	100/case
5011-300	96 well PCR tube Racks (yellow)	100/case
5011-500	96 well PCR tube Racks (orange)	100/case
5011-100	96 well PCR tube Racks (blue)	100/case
5011-200	96 well PCR tube Racks (green)	100/case

- ◆ Freezer racks with separate lids are for all major brands of screw cap micro centrifuge tubes.
- ◆ The box holds up to 0.2ml / 96 tubes.
- ◆ Has an alpha numeric labeling area.
- ◆ Have every kind of color, the convenient laboratory.
- ◆ Freezes to -90°C
- ◆ Material : Polypropylene