



81 Place Freezer Box

Cat No	specification	
5041-000	81 well Storage/Freezer Box (natural)	100/case
5041-100	81 well Storage/Freezer Box (blue)	100/case
5041-200	81 well Storage/Freezer Box (green)	100/case
5041-300	81 well Storage/Freezer Box (yellow)	100/case
5041-400	81 well Storage/Freezer Box (purple)	100/case
5041-500	81 well Storage/Freezer Box (orange)	100/case



- ◆ Freezer racks with separate lids are for all major brands of screw cap micro centrifuge tubes.
- ◆ The box holds up to 81 tubes and has an alpha numeric labeling area. Freezes to -90°C .
- ◆ Have every kind of color, the convenient laboratory.
- ◆ classification management
- ◆ Material : Polypropylene