

Multi-function box

Centrifuge tube box / Tip box

Description:

Multi function test tube can be replaced with different size test tubes or 5ml tips. The storage boxes for 15ml and 50ml centrifuge tubes can be placed. High temperature and high pressure sterilization at 121°C 15min, When the laminates are removed, they can be used as equipment sterilizing boxes or when the upper covers are removed. Can be used as small sterilize bucket or garbage bin.

- ◆ 15ml can put 5 x 5 holes and 25 test tubes. The size is W14.5 x D14.5 x H15.3 cm.
- ◆ 50ml can put 4 x 4 holes and 16 test tubes. The size is W14.5 x D14.5 x H15.3 cm.
- ◆ The 5ml tip box can be fitted with 50 tips. Can be shown in Figure tip length 14.5cm , 12.4cm or others).



SG5052-000
15ml Centrifuge tube box



SG5053-000
50ml Centrifuge tube box



SG5051-000
5ml Tip box



Cat No	specification	
SG5052-000	15ml Centrifuge tube box	20ea/cs
SG5053-000	50ml Centrifuge tube box	20ea/cs
SG5051-000	5ml Tip box	20ea/cs