


## 2.2ml 96 Well Storage Plate

- ◆ In addition to storage and screening applications, the 2.2ml Storage Plate is also
- ◆ Recommended for the growth of bacterial cultures.
- ◆ 96 square wells and 2.2ml volume.
- ◆ SBS standard footprint
- ◆ U-shaped well bottom to aid resuspension
- ◆ Material : Polypropylene



Cat No	specification	
7000-01S	96 Deep Well Square Plate, nature, Sterile	1's/bag, 50/case
7000-010	96 Deep Well Square Plate, nature	1's/bag, 50/case