

Pipette Tip

- ◆ Medical grade resin confirming to U.S.P. Class VI
- ◆ These tips fit all major pipets and offer the most versatility of any tips.
- ◆ Pipet tips with Ultrafine point feature a thin flexible section at the end that is non beveled and very flexible. This helps you dispense the last tiny amount of sample from the tip.
- ◆ Pipet tips is thin Wall · low retention pipet tips and all tips are autoclavable.
- ◆ Certified DNase & RNase Free · Pyrogen Free · Endotoxin Free · Human DNA free



0.1-10ul

Graduated at 2, 10 ul / Length 31mm



10-200ul

Graduated at 20, 200 ul / Length 50mm



100-1250ul


Length 88mm



10-200ul



10-200ul



100-1250ul



1000-5000ul

Graduated at 1.5, 2, 3, 4, 5ml / Length 14.45cm

Cat No	specification	
2011-000	0.1-10ul pipet tips · nalture	1000` s/bag, 10bag/cs
2022-000	1-200ul pipet tips · nalture	1000` s/bag, 10bag/cs
2022-300	1-200ul pipet tips · yellow	1000` s/bag, 10bag/cs
2041-000	100-1250ul pipet tips · nalture	1000` s/bag, 10bag/cs
2047-000	1000-5000ul pipet tips · nalture	250` s/bag, 10bag/cs